

PDX1 antibody



Catalog Number: 113640

Product name

PDX1 antibody

Specificity

Human; other species not tested.

Antibody description

PDX1 Rabbit Polyclonal antibody. Positive WB detected in BGC-823 cells, PC-3 cells. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000209). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



BGC-823 cells were subjected to SDS PAGE followed by western blot with Catalog No:113640(PDX1 Antibody) at dilution of 1:600