

PELI1 antibody

Catalog Number: 113650

Product name

PELI1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PELI1 Rabbit Polyclonal antibody. Positive IP detected in SH-SY5Y cells. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of PELI1 recombinant protein (Accession Number: BC011419). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:113650(PELI1 antibody) at dilution of 1:400



IP Result of anti-PELI1 (IP:Catalog No:113650, 4ug; Detection:Catalog No:113650 1:300) with SH-SY5Y cells lysate 3000ug.