

PDE1A antibody

Catalog Number: 113657

Product name

PDE1A antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

PDE1A Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 61 kDa

Preparation

This antibody was obtained by immunization of PDE1A recombinant protein (Accession Number: XM_017004295). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113657(PDE1A antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113657(PDE1A antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-PDE1A (IP:Catalog No:113657, 4ug; Detection:Catalog No:113657 1:500) with

PDE1A antibody



Catalog Number: 113657

mouse brain tissue lysate 4000ug.