

# PCF11 antibody



Catalog Number: 113703

---

## Product name

PCF11 antibody

## Immunogen

[Human PCF11 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

PCF11 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HeLa cells. Observed molecular weight by Western-blot: 173 kDa

## Preparation

This antibody was obtained by immunization of PCF11 recombinant protein (Accession Number: NM\_015885). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113703(PCF11 antibody) at dilution of 1:1000