

PCP4L1 antibody



Catalog Number: 113750

Product name

PCP4L1 antibody

Immunogen

[Human PCP4L1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCP4L1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 8-12 kDa

Preparation

This antibody was obtained by immunization of PCP4L1 recombinant protein (Accession Number: BC028905). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113750(PCP4L1 Antibody) at dilution of 1:600