

PGC1a antibody

Catalog Number: 113781

Product name

PGC1a antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PGC1a Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 91-98 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_013261). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:2000

Validations



IP Result of anti-PPARGC1A (IP:Catalog No:113781, 4ug; Detection:Catalog No:113781 1:500) with mouse skeletal muscle tissue lysate 2200ug.