

PGLYRP4 antibody

Catalog Number: 113788

Product name

PGLYRP4 antibody

Specificity

Human; other species not tested.

Antibody description

PGLYRP4 Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue, human oesophagus tissue.

Preparation

This antibody was obtained by immunization of PGLYRP4 recombinant protein (Accession Number: NM_020393). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:113788(PGLYRP4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:113788(PGLYRP4 Antibody) at dilution of 1:50 (under 10x lens)