

# PGS1 antibody

Catalog Number: 113795

## Product name

PGS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PGS1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 63 kDa

## Preparation

This antibody was obtained by immunization of PGS1 recombinant protein (Accession Number: BC025951). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113795(PGS1 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using PGS1 antibody Catalog No:113795 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).