

# PHLDA2 antibody



Catalog Number: 113824

## Product name

PHLDA2 antibody

## Specificity

human; other species not tested.

## Antibody description

PHLDA2 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, PC-3 cells. Observed molecular weight by Western-blot: 25kd

## Preparation

This antibody was obtained by immunization of PHLDA2 recombinant protein (Accession Number: NM\_003311). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113824(PHLDA2 antibody) at dilution of 1:300