

PI4K2A antibody



Catalog Number: 113876

Product name

PI4K2A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PI4K2A Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, K-562 cells. Observed molecular weight by Western-blot: 48-52kd

Preparation

This antibody was obtained by immunization of PI4K2A recombinant protein (Accession Number: NM_018425). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113876(PI4K2A antibody) at dilution of 1:800