

Piezo1 antibody



Catalog Number: 113887

Product name

Piezo1 antibody

Specificity

Human; other species not tested.

Antibody description

Piezo1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, COLO 320 cells, HeLa cells, human brain tissue, HUVEC cells. Observed molecular weight by Western-blot: 230-280 kDa

Preparation

This antibody was obtained by immunization of Piezo1 recombinant protein (Accession Number: BC008073). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113887(FAM38A antibody) at dilution of 1:400