

PKP4 antibody

Catalog Number: 113922

Product name

PKP4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PKP4 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HepG2 cells. Observed molecular weight by Western-blot: 125-140kd

Preparation

This antibody was obtained by immunization of PKP4 recombinant protein (Accession Number: NM_001005476). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113922(PKP4 antibody) at dilution of 1:300