

PLAC8 antibody

Catalog Number: 113931

Product name

PLAC8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLAC8 Rabbit Polyclonal antibody. Positive IHC detected in human tonsil tissue, human brain tissue, human placenta tissue, human spleen tissue.

Preparation

This antibody was obtained by immunization of PLAC8 recombinant protein (Accession Number: NM_016619). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human tonsil using Catalog No:113931(PLAC8 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human tonsil using Catalog No:113931(PLAC8 antibody) at dilution of 1:50 (under 40x lens)