

Placental lactogen antibody



Catalog Number: 113932

Product name

Placental lactogen antibody

Specificity

Human; other species not tested.

Antibody description

Placental lactogen Rabbit Polyclonal antibody.
Positive IHC detected in human placenta tissue.
Positive WB detected in human placenta tissue.
Observed molecular weight by Western-blot: 23-25 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001317).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113932(CSH1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human placenta using Catalog No:113932(CSH1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:113932(CSH1 antibody) at dilution of 1:50 (under 40x lens)