

# Placental lactogen antibody



Catalog Number: 113934

## Product name

Placental lactogen antibody

## Specificity

Human; other species not tested.

## Antibody description

Placental lactogen Rabbit Polyclonal antibody.  
Positive WB detected in human placenta tissue.  
Observed molecular weight by Western-blot: 22-25 kDa

## Preparation

This antibody was obtained by immunization of Placental lactogen recombinant protein (Accession Number: NM\_001317). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113934(CSH1 antibody) at dilution of 1:500