

# PLCH1 antibody

Catalog Number: 113949

## Product name

PLCH1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PLCH1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, SH-SY5Y cells. Positive IF detected in Hela. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 185kd

## Preparation

This antibody was obtained by immunization of PLCH1 recombinant protein (Accession Number: BC113950). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113949(PLCH1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:113949(PLCH1 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:113949(PLCH1

# PLCH1 antibody

Catalog Number: 113949

---

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using PLCH1 antibody Catalog No:113949 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).