

PLEKHF1 antibody



Catalog Number: 113961

Product name

PLEKHF1 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse; other species not tested.

IF: 1:10-1:100

Validations

Antibody description

PLEKHF1 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive IF detected in HepG2 cells. Positive WB detected in human skeletal muscle tissue, human kidney tissue, mouse lung tissue, mouse skeletal muscle tissue, mouse spleen tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 31kd




Preparation

This antibody was obtained by immunization of PLEKHF1 recombinant protein (Accession Number: NM_024310). Purification method: Antigen affinity purified.

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113961(PLEKHF1 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:113961(PLEKHF1 Antibody) at dilution of 1:50

Ig Type

Rabbit IgG



Applications

ELISA, WB, IF, IHC


Dilutions

Recommended Dilution:

PLEKHF1 antibody

Catalog Number: 113961

Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:113961(PLEKHF1 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells, using PLEKHF1 antibody Catalog No:113961 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).