

# PNPLA3 antibody

Catalog Number: 113981

## Product name

PNPLA3 antibody

## Specificity

Human; other species not tested.

## Antibody description

PNPLA3 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, human adipose tissue, Jurkat cells. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of PNPLA3 recombinant protein (Accession Number: BC014449). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113981(PNPLA3 antibody) at dilution of 1:100



IP Result of anti-PNPLA3 (IP:Catalog No:113981, 4ug; Detection:Catalog No:113981 1:100) with HepG2 cells lysate 1000ug.