

Podoplanin antibody



Catalog Number: 113990

Product name

Podoplanin antibody

Specificity

Human; other species not tested.

Antibody description

Podoplanin Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human placenta tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of Podoplanin recombinant protein (Accession Number: BC022812). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113990(PDPN antibody) at dilution of 1:1000