

POF1B antibody



Catalog Number: 113992

Product name

POF1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

POF1B Rabbit Polyclonal antibody. Positive WB detected in human bladder tissue, mouse bladder tissue, mouse lung tissue. Observed molecular weight by Western-blot: 69kd

Preparation

This antibody was obtained by immunization of POF1B recombinant protein (Accession Number: NM_001307940). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human bladder tissue were subjected to SDS PAGE followed by western blot with Catalog No:113992(POF1B antibody) at dilution of 1:500