

PLS1 antibody

Catalog Number: 114009

Product name

PLS1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PLS1 Rabbit Polyclonal antibody. Positive IHC detected in human colon. Positive WB detected in mouse small intestine tissue, COLO 320 cells. Observed molecular weight by Western-blot: 70kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002670). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:114009(PLS1 antibody) at dilution of 1:1200



Immunohistochemical of paraffin-embedded human colon using Catalog No:114009(PLS1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon using Catalog No:114009(PLS1 antibody) at dilution of 1:50 (under 40x lens)