

# PLSCR3 antibody



Catalog Number: 114013

## Product name

PLSCR3 antibody

## Specificity

Human; other species not tested.

## Antibody description

PLSCR3 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human skeletal muscle tissue.

## Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM\_020360). Purification method: Protein A purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114013(PLSCR3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114013(PLSCR3 antibody) at dilution of 1:50 (under 40x lens)