

PMFBP1 antibody



Catalog Number: 114023

Product name

PMFBP1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PMFBP1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HepG2 cells, MCF7 cells. Observed molecular weight by Western-blot: 119kd

Preparation

This antibody was obtained by immunization of PMFBP1 recombinant protein (Accession Number: BC028423). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114023(PMFBP1 antibody) at dilution of 1:500