

POLR1A antibody

Catalog Number: 114034

Product name

POLR1A antibody

Specificity

Human; other species not tested.

Antibody description

POLR1A Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 190-195 kDa

Preparation

This antibody was obtained by immunization of POLR1A recombinant protein (Accession Number: BC017115). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114034(RPA1 antibody) at dilution of 1:1000



IP Result of anti-RPA1 (IP:Catalog No:114034, 4ug; Detection:Catalog No:114034 1:300) with HeLa cells lysate 2400ug.