

PPAN antibody

Catalog Number: 114079

Product name

PPAN antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPAN Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive IHC detected in human lymphoma tissue. Positive WB detected in mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 55-58 kDa

Preparation

This antibody was obtained by immunization of PPAN recombinant protein (Accession Number: NM_020230). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114079(PPAN antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:114079(PPAN antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of MCF-7 cells, using Catalog No:114079 and FITC-labeled donkey anti-rabbit IgG

PPAN antibody



Catalog Number: 114079

rabbit IgG (green). Blue pseudocolor = DAPI (fluorescent DNA dye).