

PPP4C antibody

Catalog Number: 114116

Product name

PPP4C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPP4C Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, mouse kidney tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of PPP4C recombinant protein (Accession Number: NM_002720). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114116(PPP4C antibody) at dilution of 1:1000



IP Result of anti-PPP4C (IP:Catalog No:114116, 4ug; Detection:Catalog No:114116 1:1000) with mouse kidney tissue lysate 4000ug.