

PPFIA3 antibody

Catalog Number: 114120

Product name

PPFIA3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPFIA3 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, HEK-293 cells. Observed molecular weight by Western-blot: 140 kDa

Preparation

This antibody was obtained by immunization of PPFIA3 recombinant protein (Accession Number: BC021255). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114120(PPFIA3 Antibody) at dilution of 1:600



IP Result of anti-PPFIA3 (IP:Catalog No:114120, 4ug; Detection:Catalog No:114120 1:600) with mouse brain tissue lysate 3440ug.