

# PPIL2 antibody



Catalog Number: 114128

## Product name

PPIL2 antibody

## Immunogen

[Human PPIL2 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PPIL2 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, human adrenal gland tissue. Observed molecular weight by Western-blot: 59kd

## Preparation

This antibody was obtained by immunization of PPIL2 recombinant protein (Accession Number: NM\_001317996). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114128(PPIL2 antibody) at dilution of 1:400