

PPM1A antibody



Catalog Number: 114132

Product name

PPM1A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPM1A Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HeLa cells, HepG2 cells, human heart tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of PPM1A recombinant protein (Accession Number: XM_011536879). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114132(PPM1A antibody) at dilution of 1:900