

PRAM1 antibody



Catalog Number: 114164

Product name

PRAM1 antibody

Specificity

Human; other species not tested.

Antibody description

PRAM1 Rabbit Polyclonal antibody. Positive WB detected in U-937 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of PRAM1 recombinant protein (Accession Number: NM_032152). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:114164(PRAM1 Antibody) at dilution of 1:300