

# PRH1 antibody

Catalog Number: 114188

## Product name

PRH1 antibody

## Specificity

Human; other species not tested.

## Antibody description

PRH1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Transfected HEK-293 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 17 kDa

## Preparation

This antibody was obtained by immunization of PRH1 recombinant protein (Accession Number: BC064553). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200


## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114188(PRH1 Antibody) at dilution of 1:350



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114188(PRH1 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:114188(PRH1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)