

PRKG1 antibody

Catalog Number: 114202

Product name

PRKG1 antibody

Specificity

Human, Rat, Mouse; other species not tested.

Antibody description

PRKG1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in rat brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 65-78 kDa

Preparation

This antibody was obtained by immunization of PRKG1 recombinant protein (Accession Number: NM_006258). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114202(PRKG1 Antibody) at dilution of 1:600



IP Result of anti-PRKG1 (IP:Catalog No:114202, 4ug; Detection:Catalog No:114202 1:600) with mouse brain tissue lysate 4000ug.