

PSTPIP1 antibody



Catalog Number: 114255

Product name

PSTPIP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PSTPIP1 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 45-48kd

Preparation

This antibody was obtained by immunization of PSTPIP1 recombinant protein (Accession Number: NM_003978). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114255(PSTPIP1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:114255(PSTPIP1 antibody) at dilution of 1:50 (under 40x lens)