

# PTBP1 antibody



Catalog Number: 114258

## Product name

PTBP1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PTBP1 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HepG2 cells, A431 cells. Observed molecular weight by Western-blot: 55-62kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002819).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114258(PTBP1 antibody) at dilution of 1:1000



Immunofluorescent analysis of HEK-293 cells using Catalog No:114258(PTBP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)