

# PTCH1 antibody

Catalog Number: 114262

## Product name

PTCH1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PTCH1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in mouse brain tissue, A431 cells, HeLa cells. Observed molecular weight by Western-blot: 161 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000264). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114262(PTCH1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:114262(PTCH1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:114262(PTCH1 antibody) at dilution of 1:50 (under 40x lens)