

PTGER3 antibody



Catalog Number: 114267

Product name

PTGER3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PTGER3 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue.

Preparation

This antibody was obtained by immunization of PTGER3 recombinant protein (Accession Number: NM_198716). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114267(PTGER3 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114267(PTGER3 Antibody) at dilution of 1:200 (under 40x lens).