

# PTMS antibody



Catalog Number: 114308

## Product name

PTMS antibody

## Specificity

Human; other species not tested.

## Antibody description

PTMS Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive WB detected in U-937 cells. Observed molecular weight by Western-blot: 10-12kd

## Preparation

This antibody was obtained by immunization of PTMS recombinant protein (Accession Number: NM\_002824). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:114308(PTMS antibody) at dilution of 1:100



Immunofluorescent analysis of A549 cells using Catalog No:114308(PTMS Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)