

# QKI antibody

Catalog Number: 114367

## Product name

QKI antibody

## Immunogen

[Human QKI Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

QKI Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 38kd

## Preparation

This antibody was obtained by immunization of QKI recombinant protein (Accession Number: NM\_006775). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114367(QKI antibody) at dilution of 1:400



IP Result of anti-QKI (IP:Catalog No:114367, 4ug; Detection:Catalog No:114367 1:400) with K-562 cells lysate 3600ug.