

# PSMB10 antibody



Catalog Number: 114380

## Product name

PSMB10 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PSMB10 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 29kd

## Preparation

This antibody was obtained by immunization of PSMB10 recombinant protein (Accession Number: NM\_002801). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114380(PSMB10 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114380(PSMB10 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114380(PSMB10 antibody) at dilution of 1:50 (under 40x lens)