

RAB20 antibody

Catalog Number: 114419

Product name

RAB20 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB20 Rabbit Polyclonal antibody. Positive IP detected in SW 1990 cells. Positive WB detected in SW 1990 cells, HeLa cells. Observed molecular weight by Western-blot: 26-29 kDa

Preparation

This antibody was obtained by immunization of RAB20 recombinant protein (Accession Number: NM_017817). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



SW 1990 cells were subjected to SDS PAGE followed by western blot with Catalog No:114419(RAB20 antibody) at dilution of 1:1000



IP Result of anti-RAB20 (IP:Catalog No:114419, 4ug; Detection:Catalog No:114419 1:1000) with SW 1990 cells lysate 2400ug.