

RASIP1 antibody

Catalog Number: 114481

Product name

RASIP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RASIP1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, NIH/3T3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 115 kDa

Preparation

This antibody was obtained by immunization of RASIP1 recombinant protein (Accession Number: BC021860). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114481(RASIP1 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:114481(RASIP1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)