

# RABEP1 antibody

Catalog Number: 114493

## Product name

RABEP1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RABEP1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, A431 cells, human liver tissue. Observed molecular weight by Western-blot: 100 kDa

## Validations




## Preparation

This antibody was obtained by immunization of RABEP1 recombinant protein (Accession Number: NM\_004703). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114493(RABEP1 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

IP Result of anti-RABEP1 (IP:Catalog No:114493, 3ug; Detection:Catalog No:114493 1:500) with HeLa cells lysate 3800ug.

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IP, IF

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog

# RABEP1 antibody



Catalog Number: 114493

---

No:114493(RABEP1 Antibody) at dilution of 1:50

and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)