

RAC1/2/3 antibody



Catalog Number: 114504

Product name

RAC1/2/3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

RAC1/2/3 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human brain tissue. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 21 kDa

Validations



Preparation

This antibody was obtained by immunization of RAC1/2/3 recombinant protein (Accession Number: NM_005052). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114504(RAC1/2/3 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF, IP


Dilutions

Recommended Dilution:

WB: 1:500-1:5000




Immunohistochemical of paraffin-embedded human brain using Catalog No:114504(RAC1/2/3 antibody) at dilution of 1:100 (under 10x lens)



RAC1/2/3 antibody


Catalog Number: 114504

Immunohistochemical of paraffin-embedded human brain using Catalog No:114504(RAC1/2/3 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using

RAC3 antibody Catalog No:114504 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



IP Result of anti-RAC1/2/3 (IP:Catalog No:114504, 3ug; Detection:Catalog No:114504 1:500) with mouse brain tissue lysate 2640ug.