

RBMY1A1 antibody



Catalog Number: 114558

Product name

RBMY1A1 antibody

Specificity

Human; other species not tested.

Antibody description

RBMY1A1 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, NCCIT cell. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of RBMY1A1 recombinant protein (Accession Number: XM_011531490). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:114558(RBMY1A1 Antibody) at dilution of 1:300