

# RBM10 antibody



Catalog Number: 114595

## Product name

RBM10 antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RBM10 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 120 kDa

## Preparation

This antibody was obtained by immunization of RBM10 recombinant protein (Accession Number: NM\_005676). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114595(RBM10 antibody) at dilution of 1:1000



IP Result of anti-RBM10 (IP:Catalog No:114595, 4ug; Detection:Catalog No:114595 1:800) with mouse brain tissue lysate 4000ug.