

RBM17 antibody



Catalog Number: 114602

Product name

RBM17 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RBM17 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, HeLa cells, human kidney tissue, MCF7 cells. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of RBM17 recombinant protein (Accession Number: NM_001145547). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114602(RBM17 antibody) at dilution of 1:500