

RFP2 antibody

Catalog Number: 114629

General Information

Product name:

RFP2 antibody

Immunogen:

Applications:

ELISA, WB, IHC, IP

Specificity:

Human, Mouse; other species not tested.

Antibody description:

RFP2 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in Jurkat cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 47 kDa

Properties

Preparation:

This antibody was obtained by immunization of RFP2 recombinant protein (Accession Number: BC003579). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal

Ig Type:Rabbit IgG

Applications


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:500-1:5000

IHC: 1:20-1:200

Validations





Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114629(TRIM13 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:114629(TRIM13 Antibody) at dilution of 1:50

RFP2 antibody

Catalog Number: 114629



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:114629 (TRIM13 Antibody) at dilution of 1:50

IP Result of anti-TRIM13 (IP:Catalog No:114629, 4ug; Detection:Catalog No:114629 1:1000) with Jurkat cells lysate 3200ug.