

RGSL1 antibody



Catalog Number: 114651

Product name

RGSL1 antibody

Immunogen

[Human RGSL1 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

RGSL1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 125 kDa

Preparation

This antibody was obtained by immunization of RGSL1 recombinant protein (Accession Number: NM_001137669). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114651(RGSL1 Antibody) at dilution of 1:600