

RHAMM,Lgals9 antibody



Catalog Number: 114652

Product name

RHAMM,Lgals9 antibody

Specificity

mouse, rat; other species not tested.

Antibody description

RHAMM,Lgals9 Rabbit Polyclonal antibody.
Positive WB detected in NIH/3T3 cells. Observed
molecular weight by Western-blot: 32kd,70kd

Preparation

This antibody was obtained by immunization of
Peptide (Accession Number: NM_012964).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE
followed by western blot with Catalog
No:114652(RHAMM,Lgals9 antibody) at dilution
of 1:400