

RDH14 antibody



Catalog Number: 114670

Product name

RDH14 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RDH14 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, A549 cells. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of RDH14 recombinant protein (Accession Number: NM_020905). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114670(RDH14 antibody) at dilution of 1:200